



# REMINDER



## AED ELECTRODE PADS and BATTERY → EXPIRATION DATE

Dear Valued Client:

You may own a Cardiac Science **Powerheart G3 Plus Automated External Defibrillator (AED)**. The electrodes have a shelf life of two years. Each packet has an expiration date clearly marked on them. The adhesive material starts to dry out and the battery needs replaced after this allotted time. It is neither wise nor safe to use in a cardiac arrest medical emergency. Your adult and/or pediatric electrodes may have expired so please check your AED for this and for the **Green** Rescue Ready symbol on the upper right side of the AED handle signifying that the battery is also in working order. If you have a **Red** Rescue Ready symbol with a **Black X**, there is a problem with the pads and/or battery.

Please complete the order form below if you need the adult/pediatric electrode pads and non-rechargeable lithium battery. It is recommended by the AHA that you keep 2 sets of adult electrodes with each AED unit (1 attached to the AED and 1 as a spare). You may fax to the number below or scan to Scott @LifeServers.com. Prices include all sales tax and freight costs.

Thank you for your past business and for thinking of the health and well-being of your members, co-workers, or employees. Please contact us should you have any other interest in AEDs.

Please Detach Here

### CATALOG #9131-001—Adult Defibrillation Pads

\_\_\_\_\_ 2 Pair Adult Defibrillation Pads

### CATALOG #9730-002—Pediatric Defibrillation Pads

\_\_\_\_\_ 1 Pair Pediatric Defibrillation Pads

### CATALOG #9146-202—Intellisense Lithium Battery

\_\_\_\_\_ 1 Lithium Battery

COMPANY: \_\_\_\_\_

ATTN: \_\_\_\_\_

STREET: \_\_\_\_\_

CITY, STATE, ZIP: \_\_\_\_\_

PHONE #: \_\_\_\_\_

EMAIL: \_\_\_\_\_

### **Authorized Distributor for Cardiac Science and Physio-Control Defibrillators (AEDs)**

2926 Pickwick Drive | Columbus, Ohio 43221 | 614.488.6603 | Fax: 614.488.6605 | [www.lifeservers.com](http://www.lifeservers.com)