

# POWERHEART® AED G3 Accessories and Supplies

Designed to give you the most out of your Powerheart AED

## Primary Users

G3 and G3 Plus Accessories

- Health care facilities
- Sports venues
- Schools
- Recreation facilities
- Places of worship
- Any public place

G3 Pro® Accessories:

- EMS responders
- Hospitals
- Medical professionals

## Primary Benefits

**Ease of use.** Adult non-polarized electrodes (9131) can be placed in either position on the chest (upper-right or lower-left), reducing confusion and saving you valuable time during a rescue.

**Reliability.** The 9146 medical-grade Intellisense® lithium battery for the Powerheart AED G3 and G3 Plus comes with a four-year full operational guarantee.

**Flexibility for professionals.** Professional responders will appreciate additional tools for the G3 Pro including a 3-lead ECG cable, polarized electrodes, and a rechargeable battery.

**Pediatric applications.** Our pediatric electrodes reduce the energy delivered by Powerheart AED G3 models to the appropriate, pediatric dose for children up to 8 years of age or 55 lbs (25 kg).



## Electrodes for every emergency situation

Cardiac Science offers electrodes for a variety of rescue situations. We offer non-polarized electrodes that simplify AED use for inexperienced users, and pediatric defibrillator electrodes that reduce defibrillation energy to the appropriate pediatric dose. Polarized electrodes are available for use by trained medical responders.

## Batteries you can rely on

All Powerheart AEDs use a single, medical-grade Intellisense lithium battery (as opposed to the multiple, consumer-grade batteries required by some AED manufacturers). Each Intellisense battery can deliver at least nine shocks after the AED indicates the battery is low.

## Patient cables for the G3 Pro

The three-lead patient cable enables professional responders using the G3 Pro to view the patient's ECG without using defibrillation pads. Cables are available with either AHA or IEC labeling.

## Kits to support an AED rescue

Cardiac science has an array of kits — from a basic AED Ready Kit to the Total Response Rescue Kits — to support an AED rescue. The basic kit includes gloves, wipes, and CPR masks while total response kits provide emergency oxygen and a first aid kit.

## Data cables

We also offer data cables for communication between the AED and a computer making it possible to quickly upload rescue data and transfer critical information.



# POWERHEART® AED G3 Accessories and Supplies

Accessories and supplies available	
<b>ELECTRODES</b>	
9131 Adult Electrodes	These easy-to-use non-polarized adult defibrillation pads (either pad can be applied to either side) are intended for use with the G3 Plus, G3 Automatic, and G3 Semi-Automatic. Two-year shelf life.
9660 Adult Electrodes	These polarized (site-specific) adult defibrillation pads are intended for use to defibrillate and monitor with the G3 Pro. Two-year shelf life.
9730 Pediatric Electrodes	These pediatric defibrillation pads can be used with any AED G3 to reduce defibrillation energy to the appropriate pediatric dose. Two-year shelf life.
<b>BATTERIES</b>	
9146 Intellisense Lithium Battery	For use with the G3 and G3 Plus, this Intellisense lithium battery comes with a 4-year, full operational replacement guarantee.
9145 Intellisense Lithium Battery	For use with the G3 Pro, this Intellisense lithium battery comes with a 1-year from deployment or 12 hours of use (whichever comes first) full operational replacement guarantee.
9144 Rechargeable Battery	This rechargeable battery is for use with the G3 Pro. A fully charged battery delivers 60 to 100 shocks and has a usable lifetime of 2.5 years or 300 discharge cycles.
9044 Battery Charger	This charger will recharge the 9144 battery for the G3 Pro to full capacity.
<b>PATIENT CABLES</b>	
5111-001 3-Lead Cable, AHA labeling	This 3-lead ECG cable for G3 Pro allows users to see a patient's ECG on the Pro's color screen without requiring the user to place defibrillation pads. This cable comes labeled according to AHA standards. Two-year shelf life.
5111-002 3-Lead Cable, IEC labeling	This 3-lead ECG cable for G3 Pro allows users to see a patient's ECG on the Pro's color screen without requiring the user to place defibrillation pads. This cable comes labeled according to IEC standards. Two-year shelf life.
<b>RESCUE KITS</b>	
5550-005 AED Ready Kit	This basic ready kit includes nitrile gloves, razor, scissors, towel, gauze, antiseptic wipes, one-way filter mask for CPR, and carabineer attachment.
5587-001 Total Response Rescue Kit	The kit includes an emergency oxygen cylinder, oxygen tubing, oxygen regulator, personal mask with oxygen inlet, a blood-borne pathogens responder kit, blood-borne pathogens cleanup kit, and first aid kit. Includes a wire wall-mount rack so you can store the kit on the wall near your AED.
5588-001 Total Response Rescue Backpack	The sturdy backpack with adjustable straps includes an emergency oxygen cylinder, oxygen tubing, oxygen regulator, personal mask with oxygen inlet, a blood-borne pathogens responder kit, blood-borne pathogens cleanup kit, and first aid kit.
<b>DATA CABLES</b>	
162-0108-001 Infrared Communications Cable	For use with G3 Pro, this cable allows users to download rescue data to a PC running Rescuelink software.
170-2120 Serial Communications Cable	For use with all G3 and G3 Plus models, this cable allows users to download rescue data to a PC running Rescuelink software.

Some items may not be available in all markets. Please contact your local Cardiac Science representative for more information.

**Cardiac Science Corporation** • 3303 Monte Villa Parkway, Bothell, WA 98021 USA • 425.402.2000 • US toll-free 800.426.0337 • Fax: 425.402.2001 • info@cardiacscience.com  
**Orders and Customer Care (US and International)** • 425.402.2000 • US toll-free 800.426.0337 • Fax: 425.402.2001 • care@cardiacscience.com  
**Technical Support** • (US) Fax: 425.402.2022 • technicalsupport@cardiacscience.com • http://websupport.cardiacscience.com/webchat/ • (International) internationalservice@cardiacscience.com  
**Cardiac Science International A/S** • Kirke Vaerloesevej 14, DK-3500 Vaerloese, Denmark • +45.4438.0500 • Fax: +45.4438.0501 • international@cardiacscience.com  
**United Kingdom** • The Manse, 39 Northenden Road, Sale, Manchester, M33 2DH, United Kingdom • +44.161.926.0000 • uk@cardiacscience.com  
**France** • Parc de la Duranne, 565 Rue René Descartes, 13857 Aix en Provence cedex 3, France • +33.4.88.19.92.92 • france@cardiacscience.com  
**Central Europe (D, A, CH)** • Oskar-Schindler-Strasse 3, D-50769 Köln, Germany • +49.221.33734.300 • centraleurope@cardiacscience.com  
**China** • 6/F South Building, 829, Yi Shan Road, Shanghai 200233, China • +86.21.6495.9121 • china@cardiacscience.com